Description of this Thermometer

- ON/OFF Button
- Display
- Battery compartment cover
- Measuring sensor / measuring tip
- Cleaning and disinfecting area (thermometer probe only)

#### Important Safety Instructions

Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.

This device is only for use for measuring human body temperature at mouth, rectal or axillary. Do not attempt to take temperatures at other sites such as in the ear, as it may result in false readings and may lead to injury.

- Do not use this device if you think it is damaged or notice anything unusual.
- Follow instructions for use. This document provides important instructions before first use or personal hygiene.
- The minimum measurement time while the beep is heard must be maintained without exception!


Consider that different measurement locations may require continued measuring even after the beep, see section «Measuring different body temperatures and the measurement time». Do not attempt rectal measurements on persons with rectal disorder. Doing so may aggravate or worsen the disorder.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of damage to the device if supplied with cables or tubes.


- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect the device from impact and dropping.
- Do not immerse the thermometer in water deeper than 45°!
- Avoid ambient temperatures above 30 °C. NEVER put this device where the product was purchased, or your local MicroLife service. You may contact your local MicroLife service through our website: www.microlife.com/sr/support
- Use only the commercial disinfectants listed in the section «Cleaning and Disinfecting» to clean the device to avoid damage to the device.

We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your local MicroLife-Service to arrange the test.

**WARNING:** The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only. Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Read the instructions carefully before using this device.

Tip BF applied part

**Turning on the Thermometer**
To turn on the thermometer, press the ON/OFF button 1, a short beep signals «thermometer ON». A display test is performed. All segments should be displayed.

The last measurement reading will be shown on the display 2 automatically for 2 seconds with the «m» icon. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «C» appear at the display field 2. The thermometer is now ready for use.

#### Function Test

Correct function test of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccurate), this is indicated by «ERR» on the display and a measurement becomes impossible. In this case, the thermometer must be replaced. Using the Thermometer

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «C» symbol flashes. If the beep is heard 10 times and the «C» is no longer flashing, the predictive end-temperature has been determined and the thermometer can be read off. 35 short beeps will sound when the temperature is higher than 37.5 °C in order to alert the patient that he/she may have fever. Reference: Oral temperature.

To achieve comparable results allow a 1 minute interval between measurements. Replace the battery life, turn off the thermometer by briefly pressing the ON/OFF button 1. Otherwise the thermometer will automatically turn off after about 10 minutes.

#### Measuring methods / Normal body temperature

- In the armpit (axillary) / 34.7 - 37.3 °C**
Wipe the armpit with a dry towel. Place the measuring sensor 4 under the arm into the center of the armpit so the tip is touching the skin and position the patient's arm next to the patient's body. This ensures that the arm does not affect the reading. Because the axillary takes more time to reach its stable temperature wait at least **5 minutes**, regardless of the beep sound.

- In the mouth (oral) / 35.5 - 37.5 °C**
Do not eat or drink anything hot or cold 10 minutes before the measurement. The measuring sensor should remain closed up to 12 minutes before starting a reading.

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor 4 must be in good contact with the tissue. Close your mouth and breathe the evening air through the nose to prevent the measurement from being influenced by inhaled/exhaled air. If this is not possible due to blocked airways, another method for measuring should be used.

#### Approx. measuring time: 10 seconds!

Carefully insert the measuring sensor 4 of the thermometer 2 to 3 cm into the anal aperture. The use of a probe cover and the use of a lubricant is recommended. If you are unsure of this measurement method, you should consult a professional for guidance.

#### Approx.measuring time: 10 seconds!

**Cleaning and Disinfecting**
For disinfection in home use environment, use a 70% Isopropyl alcohol swab, or a cotton tissue moistened with 70% Isopropyl alcohol to wipe surface pollutants off the thermometer probe (note: consider the application and safety instruction of the disinfectant manufacturer). Always start wiping from the end of the thermometer probe (approx. at the middle of the thermometer towards the thermometer tip). Afterwards the entire thermometer probe (see number 3 in the drawing) should be immersed in 70% Isopropyl alcohol for at least 5 minutes (max. 24 hours). Afterward, let the disinfectant dry off for 1 minute before next use. Avoid immersing or wiping the display to protect it from fading.

The thermometer is not intended for professional use. **Battery Replacement**
When the «▼» symbol (upside-down triangle) appears in the display, the battery is at low and needs replacing. To replace the battery remove the battery compartment cover 3 from the thermometer. Insert the new battery with the «+» at the top. Make sure you have a battery of the same type to hand. Batteries can be purchased at any electrical store.

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